

THE TRUE SPECIALIST IN TMS THERAPY

Global No.1 Electroceutical Company



FDA가 인정한 차별화된 기술력을 통해 글로벌 스탠다드로 앞장서겠습니다!

전 세계 환자들의 건강한 삶에 기여할 수 있는 핵심기술 개발을 위해 리메드는 지속적으로 R&D에 투자하고, 기업역량을 강화하고 있습니다. 이를 바탕으로 국내는 물론 해외에서도 그 전문성을 인정받고 있으며, 검증된 가치에 의한 신뢰를 바탕으로 더욱 다양한 솔루션 제공을 약속드립니다.



2003년 창립된 리메드는 지난 20여년간 우울증의 근본적인 해결을 위한 뇌 자극 제품의 연구개발에 힘써왔으며, 그 결과 최첨단 TMS 장비를 성공적으로 도입했습니다.

리메드 ALTMS는 국내 최초로 미국 FDA와 식약처에서 우울증으로 승인을 받았으며, 여기에 멈추지 않고 끊임없는 노력을 통해 ALTMS24, BrainStim, 그리고 BrainStim-Home(wearable TMS)에 이르는 다양한 제품을 개발하였습니다.

앞으로도 뇌 자극 분야의 선두주자로서 리메드만의 가치있는 도전과 창조 정신으로 혁신적 기술의 성장을 이루겠습니다.

전자약 분야 전문기업 리메드가 치료의 미래를 이끌어 갑니다!





왜 TMS치료를 선택해야 하는가?

TMS(Transcranial Magnetic Stimulation)는 안전하고 비침습적인 치료로서 뇌의 뉴런을 자극하기 위해 빠르게 변화하는 자기장을 이용합니다.

이는 약물치료로 효과를 보지 못했거나, 약물을 투여하지 못하는 분들에게 효과적인 치료로 제안될 수 있으며, 더 나아가 폭 넓은 치료를 위한 새로운 옵션이 될 것입니다.

TMS를 이용한 다양한 논문들이 발표되었고, 현재도 많은 연구들이 활발히 진행되고 있습니다!

- 정신건강의학과- 우울, 강박, 불안, PTSD, 인지 저하 등
- 신경과- 두통, 이명, 수면장애, 신경병성 통증, 섬유근육통 등
- 재활의학과- 운동 기능 저하, 파킨슨병, 삼킴 장애, 구음 장애 등

20년 간의 축적된 가치로 완성한 리메드 TMS 책자를 만나보세요

국내 최초의 전자약 기업 리메드가 창립 20주년을 맞아 TMS의 A to Z를 총체적으로 담아낸 가이드북을 발간했습니다.

- 비침습적 뇌 자극, TMS의 역사 및 기전
- 국내 유수 교수진들의 질환 별 근거기반 프로토콜
- 개원가에서의 다양한 활용 방안 및 필드의 노하우

【구매 문의】 1588-7395

• 과별 최신 논문 발췌

TROKE Topics i	in Stroke Rehabilitation	Prontiers Onional Research in Psychiatry as Estimated as an array as Estimated as a control of the and array
Effects of magnetic electrical recovery i prelimina	I when IRENTING before Journal wave provide a confidence and the confidence and the configuration of the configura	Three Weeks of rTMS Treatment Maintains Clinical Improvement But Not Electrophysiological Changes in Patients With Depression: A 6-Week Follow-Up Pilot Study
Demirdal, Tüla Fazıl Mustafa (y Kurt Incesú, Ozgur Tosun, Hikmet Kocylgit, Galip Akhan & Selal	Kyung Mook Chol ^{13,4} , Soo-Hee Chol ⁴ , Sang Min Lee ^{3,4} , Kuk-In Jang ^{13,4} and Jeong-Ho Chae ^{2,4,4}
Somato Mone Brouter Somato	Sensory & Motor Research	A Machine Learning Approach to Predict rTMS Therapy Response in Major Depressive Disorder Molummal Shams
ls high-fre stimulatio superior t prefrontal	quency repetitive transcranial magnetic n of the left primary motor cortex o the stimulation of the left dorsolateral cortex in fibromyalgia syndrome?	ISSN (Print) 1226-1750 ISSN (Daling 2233-8656 Imp://doi.org/10.42333MAG.2020.25.4.691
Elif Umay Altas, Ayhan Askin, Lütfullah Beşiroğlu & Aliye Tosun		Effects of 10 Hz High Frequency Repetitive Transcranial Magnetic Stimulation Combined with Somatosensory Training for Recovery of Upper Limb Motor and Hand Function in Chronic Stroke Patients
Cal Psychopharmacology and Neuroscience 2019;17(3):409-414 Copyright© 2019, Korean College of Neuropsychopharmacology		Jung-Woo Jeong ¹ , Man-Seek Han ^{2,*} , and Bo-Kyoung Song ^{1,*}
e Efficacy of Miniaturized Repetitive Transcranial Magnetic imulation in Patients with Depression		ISSN (Prop. 1226-1790) ISSN (Prop. 1226-1790) Journal of Magnetics 26(4), 467-474 (2021) Issuer (100-0220)/MAG 2021 2.66-4867
ngmin Lee ^{1,2} , Kuk-In Jang	² , Sejin Yoon ¹ , Jeong-Ho Chae ^{1,2,4}	
nal of Magnetics 26(1), 121-128	ISSN (Piting) 1226-1730 ISSN (Value) 2234-656 https://doi.org/10.4283/JMAG.2021.26.1.121	Effects of Neuromuscular Electrical Stimulation with 5 Hz High Frequency Repetitive Transcranial Magnetic Stimulation on Cerebral Activity in Chronic Stroke Patients (Randomized Controlled Trial)
Comparison of th	a Effact of Danatitiva Transcranial Magnetic Stimulation	Jung-Woo Jeong ¹ , Sung-Ryong Ma ^{1*} , and Bo-Kyoung Song ^{1*}
by Frequency	on Upper Limb Function in Acute Stroke Patients: A Randomized Controlled Trial	DBria
Hyun Gyu Cha* Department of Physical Therapy, loongbu University, Gramsan 32713, Republic of Korea		ISSN (9444) 1236-1759 ISSN (9444) 1236-1759 ISSN (9444) 1231-1459 Imperides and Magnetics 27(4), 452-458 (2022)
I Nacomital J. 2027 Mar; M(f) ad L(bid arg/10.2021 K ad N 1070-1072 Mg/ha. 2021 K ad N 1070-1072 Mg/ha. 2021 K ad ecial Review Citat Keyaman	Safety Review for Clinical Application of Repetitive Transcranial Magnetic Stimulation	A Study of Different Frequencies of Repetitive Transcranial Magnetic Stimulation on Upper Extremity Muscle Activity and Hand Function in Chronic Stroke Patients Jung-Woo Jeong', Bockyoung Song'', and Byung-II Yang'' Bockward (2000) 2000 Bock and a Edemondation (2000) 2000 Decision for soundation of Edemondation
~~~~	Won-Scok Kim, Nam-Jong Palk	Didili Stillituidiloii ELSEVIER journal homepage: http://www.journals.elsavier.com/brain-stimulation
frontiers   Frontiers in S	vitems Neuroscience Tris. Original Reasesh Processo 10 March 202 ed. 33.3589/ega 2633.20077	The effects of repetitive transcranial magnetic stimulation on body weight and food consumption in obese adults: A randomized controlled study Settors (m ² , lawb Quarg ^{bb} , Tac-Hone Kim ² , Source Hone Im ² Vanetook Kim ²
Dreck for updates	Machine learning approaches and non-linear processing of	Torong All Eur., ', Yun Ahl Le'           Boolnet 27 fiburn, 2007
airr a Karnacajan, state Health Sciences University, d States	extracted components in frontal region to predict rTMS	ORIGINALARTICLE WILEY
ND BY ss I. Dimitriadis, timent of Clinical Psychology and obiology, Faculty of Psychology, raily of Bancelona, Spain w F. Leuchter, rely of California. Los Angeles	treatment response in major depressive disorder	The effects of high-frequency repetitive transcranial magnetic stimulation on resting-state functional connectivity in obese adults
		Se-Hong Kim MD ¹¹ 9   Bo-Yong Park MSc ²³¹ 0   Kyoungseob Byeon MSc ²³   Hyunjin Park PhD ²⁴   Youngkook Kim MD ⁵   Young-Mi Eun MD ¹   Ju-Hye Chung MD ⁶

### FDA KFDA.

06

# ALTMS

다양한 질환에 특화된 코일 적용 및 사용자의 간편성을 극대화하여 합리적인 가치를 선보입니다!

국내 전자약 기업 최초 FDA 승인

# I ココ 6

P

### 특징 및 장점

- TBS를 포함한 다양한 자극모드 기본옵션 (Single, Repetitive, TBS, Pair)
- 마그네틱 샷 카운트 제한 및 소모품이 없는 최적의 효율성
- 24시간 연속 사용이 가능한 특허 받은 오일 냉각 시스템
- 터치 스크린 형식의 사용자 친화적(user-friendly) 조작 방식
- 고급스러운 외관으로 인테리어 효과 증대
- 국내 제조사로서 신속한 AS 및 응대 가능

### 기술 사양

- 전류 파형 : 대칭 이상파 (Symmetric biphasic wave)
- 펄스폭(pulse width) : 430µsec
- 펄스 주파수 : 0.1~100pps
- 자계생성부 출력 : 최대 2.5 T-peak
- 타이머 설정 범위 : 1~60 분

### 적응증

- Major depressive disorder(MDD) Obsessive compulsive disorder(OCD) Smoking cessation Anxiety disorder
- PTSD
- Addiction
- Obesity
- Sleep Disorder

Basic











- Parkinson's disease
- Stroke
- Aphasia
- Cognitive impairment
- Dysphagia
- Neuropathic pain
- Fibromyalgia



Option



# BrainStim

개별 환자 관리가 가능한 고급형 모델로 다양한 확장기능을 통해 사용성이 높습니다!

08

### 특징 및 장점

- 환자 DB 저장이 가능해져 환자 관리 용이
- 한글 지원 및 직관적인 이미지를 통한 사용성 증대
- 질환에 따라 선택적으로 적용가능한 다양한 각도의 코일 (120°, 150°, 180°)
- 뉴로네비게이션 및 로봇암과 결합하여 확장기능 사용 가능
- 마그네틱 샷 카운트 제한 및 소모품이 없는 최적의 효율성
- 24시간 연속 사용이 가능한 특허 받은 오일 냉각 시스템
- 터치 스크린 형식의 사용자 친화적(user-friendly) 조작 방식
- 고급스러운 외관으로 인테리어 효과 증대
- 국내 제조사로서 신속한 AS 및 응대 가능

### 기술 사양

- 전류 파형 : 대칭 이상파 (Symmetric biphasic wave)
- 펄스폭(pulse width) : 360µsec
- 펄스 주파수 : 0.1~100pps
- 자계생성부 출력 : 최대 3 T-peak
- 자극 시간(On Time) : 0.1 ~ 60sec
- 휴지 시간(Off Time) : 1 ~180sec

### 적응증

• PTSD

BrainStim

- Major depressive disorder(MDD) Obsessive compulsive disorder(OCD) Smoking cessation Anxiety disorder Addiction Obesity
- Sleep Disorder



# • TBS를 포함한 다양한 자극모드의 기본옵션 (Single, Repetitive, TBS, Pair)









### BrainStim

Tinnitus

Parkinson's disease

Stroke

Aphasia

Cognitive impairment

Dysphagia

Neuropathic pain

Fibromyalgia



Option

# BrainStim-Home

컴팩트한 사이즈의 웨어러블 기기로 치료의 다양성을 높여드립니다!

Rextour

# 특징 및 장점

- 강력하고 컴팩트한 헬멧형 TMS로 최적화된 치료 기대
- 컴팩트한 크기로 공간 활용 증대
- 간편한 터치 형식으로 5가지 치료모드 선택
- 병원에 자주 내원할 수 없는 환자에게 재택치료의 편리성 제공

### 기술 사양

- 자기 에너지 : 2.3 테슬라 ± (20%)
- 반복률 : 1~20Hz
- 타이머 설정 : 60분
- 설정 모드 : 1~5
- 무게 : 헬멧 1kg, 기계 본체 : 4kg
- 전원 공급 장치 입력 전력 : 220-240V~, 50/60Hz
- 전력 소비량 : 300VA
- 치수 (W*H*D) : 231*173.5*323mm

### 적응증

Stress

Panic

- Depression
- Learning ability

- Epileptic Seizure Disorder
- cognitive function
- Tinnitus
- Stroke



본 제품은 의료기기 품목신고 및 허가 이전의 제품입니다

Anxiety

• Bipolar disorder

• Alzheimer's disease

Fatigue

Hallucination

• Dementia

Pressure



#### www.remed.kr

Sales Division / R&D Center 1 경기도 성남시 수정구 위례서일로 1길 21-7 (창곡동) 산해빌딩 (주)리메드 4, 5층 Tel. 1588-7395 Fax. 031-606-5757

#### Factory 1

대전광역시 유성구 테크노2로 187 (용산동) 미건테크노월드 301~303호 Tel. 042-934-5560 Fax. 042-934-5562

**R&D Center 2 / Factory 2** 충청북도 청주시 흥덕구 오송읍 오송생명2로 84 Tel. 043-231-8921 Fax. 043-231-8923